


Product Datasheet

Mouse sIL-6R (Soluble Interleukin 6 Receptor) Quick ELISA Kit (orb2567817)

Description	Mouse sIL-6R (Soluble Interleukin 6 Receptor) Quick ELISA Kit
Reactivity	Mouse
Range	12.5-800pg/ml
Target	sIL-6R
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Uniprot ID	P22272
Sensitivity	7.5pg/ml
Expiration Date	Please enquire.



Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)