

Product Datasheet

CRMP2 (phospho-S522) peptide (orb256789)

Description	The peptide is used to anti CRMP2 (phospho-S522) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	DPYSL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256499 (CRMP2 (phospho-S522) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O08553, P47942, Q16555
Entrez	1808, 25416, 12934
Expiration Date	12 months from date of receipt.