

Product Datasheet

EGFR (phospho-T693) peptide (orb256797)

Description	The peptide is used to anti EGFR (phospho-T693) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	EGFR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256507 (EGFR (phospho-T693) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q01279 , P00533
Entrez	1956 , 13649
Expiration Date	12 months from date of receipt.