


## Product Datasheet

# Human ERK (ELK-Related Tyrosine Kinase) Quick ELISA Kit (orb2568086)

<b>Description</b>	Human ERK (ELK-Related Tyrosine Kinase) Quick ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625-40ng/ml
<b>Target</b>	ERK
<b>Storage</b>	2 - 8 °C for 12 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
<b>Uniprot ID</b>	<b>P29323</b>
<b>Sensitivity</b>	0.375ng/ml
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)