

Product Datasheet

Estrogen Receptor beta 2 peptide (orb256817)

Description	Estrogen Receptor beta 2 Peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ESR2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256527 (ESR2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O08537 , Q92731 , Q62986
Entrez	13983 , 25149 , 2100
Expiration Date	12 months from date of receipt.