

Product Datasheet

FRS2 (phospho-Y436) peptide (orb256831)

Description	The peptide is used to anti FRS2 (phospho-Y436) antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	FRS2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256541 (FRS2 (phospho-Y436) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8WU20 , Q8C180
Entrez	327826 , 10818
Expiration Date	12 months from date of receipt.