

Product Datasheet

GATA1 (phospho-S310) peptide (orb256842)

Description	The peptide is used to anti GATA1 (phospho-S310) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep
Tested Applications	Blocking
Target	GATA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256552 (GATA1 (phospho-S310) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P17679 , P43429 , P15976
Entrez	14460 , 100911167 , 2623 , 25172
Expiration Date	12 months from date of receipt.