

Product Datasheet

GATA3 (phospho-S308) peptide (orb256844)

Description	The peptide is used to anti GATA3 (phospho-S308) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Zebrafish
Tested Applications	Blocking
Target	GATA3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256554 (GATA3 (phospho-S308) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P23772 , P23771
Entrez	2625 , 14462
Expiration Date	12 months from date of receipt.