

Product Datasheet

Connexin 43 (phospho-S368) peptide (orb256848)

Description	The peptide is used to anti Connexin 43 (phospho-S368) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GJA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256558 (Connexin 43 (phospho-S368) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P08050 , P17302 , P23242
Entrez	24392 , 14609 , 2697
Expiration Date	12 months from date of receipt.