

## Product Datasheet

### GLUR2 (phospho-S880) peptide (orb256862)

<b>Description</b>	The peptide is used to anti GLUR2 (phospho-S880) antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Primate, Rat, Zebrafish
<b>Tested Applications</b>	Blocking
<b>Target</b>	GRIA2
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb256572 (GLUR2 (phospho-S880) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P42262</a> , <a href="#">P23819</a> , <a href="#">P19491</a>
<b>Entrez</b>	<a href="#">2891</a> , <a href="#">14800</a> , <a href="#">29627</a>
<b>Expiration Date</b>	12 months from date of receipt.