

Product Datasheet

GLURD2 peptide (orb256863)

Description	The peptide is used to anti GLURD2 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	GRID2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256573 (GLURD2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61625 , Q63226 , O43424
Entrez	79220 , 14804 , 2895
Expiration Date	12 months from date of receipt.