

Product Datasheet

Human Lipocalin-2(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit (orb2568735)

| | |
|------------------------|---|
| Description | Human Lipocalin-2(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit |
| Reactivity | Human |
| Range | 31.25-2000pg/ml |
| Target | Lipocalin-2 |
| Storage | 2 - 8 °C for 12 months |
| Note | For research use only |
| Uniprot ID | P80188 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)