

Product Datasheet

Histone Deacetylase 6 (phospho-S22) peptide (orb256877)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | The peptide is used to anti Histone Deacetylase 6 (phospho-S22) antibody reactivity. |
| Reactivity | Human, Primate |
| Tested Applications | Blocking |
| Target | HDAC6 |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb256587 (Histone Deacetylase 6 (phospho-S22) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | Q9UBN7 , Q9Z2V5 |
| Entrez | 10013 , 15185 |
| Expiration Date | 12 months from date of receipt. |