

## Product Datasheet

# Mouse IgG2b $\kappa$ Isotype Control, + Antibody (orb2568791)

<b>Description</b>	Mouse IgG2b $\kappa$ Isotype Control is a mouse monoclonal antibody ideal for use as a non- reactive, isot...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	In vivo
<b>Immunogen</b>	recombinant N-Protein od SARS-CoV-2 Virus
<b>Preservatives</b>	Clear Liquid, PBS, pH 7.4, 0.2 $\mu$ m sterile filtered
<b>Form/Appearance</b>	Clear Liquid, PBS, pH 7.4, 0.2 $\mu$ m sterile filtered
<b>Concentration</b>	$\geq$ 0.5 mg/ mL
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b $\kappa$
<b>Clonality</b>	Monoclonal
<b>Purity</b>	$\geq$ 90% (CGE, reducing conditions). $\leq$ 10 % aggregates (analytical SEC). Endotoxin level $\leq$ 1 EU/mg
<b>Expiration Date</b>	12 months from date of receipt.