

Product Datasheet

SHIP1 (phospho-Y1021) peptide (orb256918)

Description	The peptide is used to anti SHIP1 (phospho-Y1021) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	INPP5D
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256628 (SHIP1 (phospho-Y1021) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q92835, Q9ES52, P97573
Entrez	54259, 3635, 16331
Expiration Date	12 months from date of receipt.