

Product Datasheet

JUNB (phospho-S259) peptide (orb256931)

| | |
|----------------------------|---|
| Description | The peptide is used to anti JUNB (phospho-S259) antibody reactivity. |
| Reactivity | Human, Mouse, Rat, Sheep |
| Tested Applications | Blocking |
| Target | JUNB |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb256641 (JUNB (phospho-S259) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | P24898 , P09450 , P17275 |
| Entrez | 16477 , 24517 , 3726 |
| Expiration Date | 12 months from date of receipt. |