

Product Datasheet

MAF1 peptide (orb256954)

Description	The peptide is used to anti MAF1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MAF1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256664 (MAF1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9H063 , Q9D0U6 , Q5XIH0
Entrez	315093 , 68877 , 84232
Expiration Date	12 months from date of receipt.