

Product Datasheet

MEKK1 (phospho-T1402) peptide (orb256959)

Description	The peptide is used to anti MEKK1 (phospho-T1402) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	MAP3K1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256669 (MEKK1 (phospho-T1402) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13233 , Q62925 , P53349
Entrez	116667 , 4214
Expiration Date	12 months from date of receipt.