

Product Datasheet

MRPS18A peptide (orb256982)

Description	The peptide is used to anti MRPS18A antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MRPS18A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256692 (MRPS18A antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9NVS2 , Q99N85
Entrez	68565 , 55168
Expiration Date	12 months from date of receipt.