

## Product Datasheet

### p47 phox (phospho-S328) peptide (orb256997)

<b>Description</b>	The peptide is used to anti p47 phox (phospho-S328) antibody reactivity.
<b>Reactivity</b>	Human, Mouse, Porcine, Primate, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	NCF1
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb256707 (p47 phox (phospho-S328) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<b>P14598, Q09014</b>
<b>Entrez</b>	<b>17969, 653361</b>
<b>Expiration Date</b>	12 months from date of receipt.