

Product Datasheet

NF-kappaB p105/p50 (phospho-S337) peptide (orb257009)

Description	The peptide is used to anti NF-kappaB p105/p50 (phospho-S337) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	NFKB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256719 (NF-kappaB p105/p50 (phospho-S337) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P19838 , P25799 , Q63369
Entrez	18033 , 4790
Expiration Date	12 months from date of receipt.