

Product Datasheet

Cystatin C antibody (FITC) (orb257054)

Description	Cystatin C antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide (please contact us for details)
Target	Cystatin C (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	15 kDa
Uniprot ID	P01034
Entrez	1471
Dilution Range	IF/ICC:1:50-200 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.