
Product Datasheet

SEMA3A antibody (FITC) (orb257060)

Description

SEMA3A antibody (FITC)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

FITC

Tested Applications

ICC, IF

Immunogen

Synthetic peptide (please contact us for details)

Target

SEMA3A (FITC)

Form/Appearance

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

86 kDa

Uniprot ID[Q14563](#), [O08665](#), [Q63548](#)**NCBI**[NP_006071.11](#), [NP_006071.1](#),
[NM_006080.2](#), [NM_006080.21](#)**Entrez**[10371](#)**Dilution Range**

IF/ICC:1:50-200 (based on 0.5 mg/ml)

Expiration Date

12 months from date of receipt.