


Product Datasheet

Human BMPR-1A / ALK-3 Protein (orb257183)

Description	Human BMPR-1A, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Arg 152 (A...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	BMPR-1A
Preservatives	50 mM Tris, 100 mM Glycine, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.3 kDa. The protein migrates as 56 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	CAA80257
Purity	98%
Source	Human BMPR-1A, Fc Tag (orb257183) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Arg 152 (Accession # P36894-1).
MW	40.3 kDa
Uniprot ID	P36894
NCBI	CAA80257
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human BMPR-1A prote...