

Product Datasheet

Human Angiopoietin-2 / ANGPT2 (275-496) Protein (orb257185)

Description

Human Angiopoietin-2 (275-496), His Tag is



contains AA L...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Angiopoietin-2

Preservatives

20 mM MOPS, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

N-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 26.3 kDa. The protein migrates as 30-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAI26201.1

Purity

95%

Source


Human Angiopoietin-2 (275-496), His Tag (orb257185) is expressed from human 293 cells (HEK293). It contains AA Lys 275 - Phe 496 (Accession # AAI26201.1).

MW

26.3 kDa

Uniprot ID
[O15123](#)
NCBI
[AAI26201.1](#)
Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human Angiopoietin-
...