


Product Datasheet

Human Angiopoietin-like 3 / ANGPTL3 Protein (orb257186)

Description

Human Angiopoietin-like 3, His Tag is expressed


 SDS-PAGE
analysis of
Human
Angiopoietin-
...

Ser 17...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Angiopoietin-like 3

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.6 kDa. The protein migrates as 28-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_055310

Purity

95%

Source

Human Angiopoietin-like 3, His Tag (orb257186) is expressed from human 293 cells (HEK293). It contains AA Ser 17 - Pro 220 (Accession # NP_055310).

MW

24.6 kDa

Uniprot ID
[Q9Y5C1](#)
NCBI
[NP_055310](#)
Expiration Date

6 months from date of receipt.