


Product Datasheet

Human Annexin A5 Protein (orb257188)

Description	Human ANXA5, Tag Free is expressed from E. coli cells. It contains AA Met 1 - Asp 320 (Accession # P...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Annexin A5
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 35.9 kDa. The protein migrates as 33-35 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH01429
Purity	98%
Source	Human ANXA5, Tag Free (orb257188) is expressed from E. coli cells. It contains AA Met 1 - Asp 320 (Accession # P08758-1).
MW	35.9 kDa
Uniprot ID	P08758
NCBI	AAH01429
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Annexin A5
pr...