



Product Datasheet

Human FGF basic Protein (orb257225)

Description	Human FGF basic, premium grade is expressed from E. coli cells. It contains AA Pro 143 - Ser 288 (Ac...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	FGF basic
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 24 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 16.5 kDa. The protein migrates as 17 kDa±3 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).
Purity	95%
Source	Human FGF basic, premium grade (orb257225) is expressed from E. coli cells. It contains AA Pro 143 - Ser 288 (Accession # P09038-4).
MW	16.5 kDa
Uniprot ID	P09038
NCBI	AAI66646.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human FGF basic pro...



Immobilized Human FGF basic, premium gra...