



Product Datasheet

Human BMP-2 Protein (orb257227)

Description	Human BMP-2 Protein, premium grade is expressed from E. coli cells. It contains AA Ala 284 - Arg 396...
Reactivity	Human
Endotoxins	0.2 EU per µg
Conjugation	Unconjugated
Target	BMP-2
Preservatives	100 mM Acetic Acid, pH2.8
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 12.8 kDa. the protein migrates as 14 kDa under reducing (R) condition, and 28 kDa under non-reducing (NR) condition (SDS-PAGE).
Purity	95%
Source	Human BMP-2 Protein, premium grade (orb257227) is expressed from E. coli cells. It contains AA Ala 284 - Arg 396 (Accession # NP_001191).
MW	12.8 kDa
Uniprot ID	P12643
NCBI	NP_001191
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human BMP-2 protein



Immobilized Human BMP-2 Protein, premium...