


Product Datasheet

Human CD3 delta Protein (orb257320)

Description	Human CD3 delta Protein
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD3 delta
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 11.4 kDa. The protein migrates as 18-28 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_000723.1
Purity	90%
Source	Human CD3 delta, His Tag (orb257320) is expressed from human 293 cells (HEK293). It contains AA Phe 22 - Ala 105 (Accession # P04234-1).
MW	11.4 kDa
Uniprot ID	P04234
NCBI	NP_000723.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human CD3 delta pro...