


Product Datasheet

Human Cystatin S / CST4 Protein (orb257368)

Description	Human Cystatin S, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 21 - Ala 14...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Cystatin S
Preservatives	25 mM MES, 150 mM NaCl, pH6.0
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.0 kDa. The protein migrates as 15 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH65714
Purity	95%
Source	Human Cystatin S, His Tag (orb257368) is expressed from human 293 cells (HEK293). It contains AA Ser 21 - Ala 141 (Accession # AAH65714).
MW	15.0 kDa
Uniprot ID	P01036
NCBI	AAH65714
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Cystatin S pr...