


Product Datasheet

Human Ephrin-A4 / EFNA4 Protein (orb257409)

Description	Human Ephrin-A4, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Gly 171 ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Ephrin-A4
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 43.0 kDa. The protein migrates as 45-50 kDa and 66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAI07484
Purity	95%
Source	Human Ephrin-A4, Fc Tag (orb257409) is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Gly 171 (Accession # AAI07484).
MW	43.0 kDa
Uniprot ID	P52798
NCBI	AAI07484
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Ephrin-A4 pro...