

Product Datasheet

SLITRK3 antibody (orb25743)

Description	SLITRK3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human SLIT and NTRK-like protein 3 protein (746-977AA)
Target	SLITRK3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O94933
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.