


Product Datasheet

Human ErbB4 / Her4 Protein (orb257443)

Description	Human ErbB4, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 26 - Pro 651 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ErbB4
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 96.0 kDa. The protein migrates as 90-106 kDa and 116-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_005226.1
Purity	90%
Source	Human ErbB4, Fc Tag (orb257443) is expressed from human 293 cells (HEK293). It contains AA Gln 26 - Pro 651 (Accession # Q15303-1).
MW	96.0 kDa
Uniprot ID	Q15303
NCBI	NP_005226.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human ErbB4 protein