


Product Datasheet

Human Fetuin A / FETUA Protein (orb257467)

Description	Human Fetuin A, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 19 - Val 367 ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Fetuin A
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 38.2 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001613.2
Purity	98%
Source	Human Fetuin A, His Tag (orb257467) is expressed from human 293 cells (HEK293). It contains AA Ala 19 - Val 367 (Accession # P02765-1).
MW	38.2 kDa
Uniprot ID	P02765
NCBI	NP_001613.2
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Fetuin A prot...