


Product Datasheet

Human Fetuin B Protein (orb257468)

Description	Human Fetuin B, His Tag is expressed from human 293 cells (HEK293). It contains AA Cys 16 - Pro 382 ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Fetuin B
Preservatives	50 mM Tris, 150 mM NaCl, pH8.0
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 42.4 kDa. The protein migrates as 53-58 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH69670
Purity	95%
Source	Human Fetuin B, His Tag (orb257468) is expressed from human 293 cells (HEK293). It contains AA Cys 16 - Pro 382 (Accession # AAH69670).
MW	42.4 kDa
Uniprot ID	Q9UGM5
NCBI	AAH69670
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Fetuin B prot...