


## Product Datasheet

### Human FGF R4 / CD334 Protein (orb257472)

<b>Description</b>	Human FGF R4, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 22 - Asp 369 (A...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	FGF R4
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-6×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.3 kDa. The protein migrates as 50-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_002002.3
<b>Purity</b>	95%
<b>Source</b>	Human FGF R4, His Tag (orb257472) is expressed from human 293 cells (HEK293). It contains AA Leu 22 - Asp 369 (Accession # P22455-1).
<b>MW</b>	39.3 kDa
<b>Uniprot ID</b>	<b>P22455</b>
<b>NCBI</b>	<b>NP_002002.3</b>
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human FGF R4 protei...