


Product Datasheet

Human Flt-3 Ligand Protein (orb257476)

Description	Human Flt-3 Ligand, premium grade is expressed from human 293 cells (HEK293). It contains AA Thr 27 ...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	Flt-3 Ligand
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 18.0 kDa. The protein migrates as 24 kDa, 27 kDa±3 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human Flt-3 Ligand, premium grade (orb257476) is expressed from human 293 cells (HEK293). It contains AA Thr 27 - Pro 185 (Accession # AAI26294).
MW	18.0 kDa
Uniprot ID	P49771
NCBI	AAI26294
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Flt-3 Ligand ...