
Product Datasheet


Human Follistatin-like 1 / FSTL1 Protein (orb257484)

Description

Human Follistatin-like 1, His Tag is expressed from



Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Follistatin-like 1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 33.8 kDa. The protein migrates as 40-55 kDa under reducing (R) condition (SDS-PAGE) due to different glycosylation.
Purity	95%
Source	Human Follistatin-like 1, His Tag (orb257484) is expressed from human 293 cells (HEK293). It contains AA Glu 21 - Ile 308 (Accession # AAH00055).
MW	33.8 kDa
Uniprot ID	Q12841
NCBI	AAH00055
Expiration Date	6 months from date of receipt.


 SDS-PAGE analysis of Human Follistatin-
I...