


Product Datasheet

Human GAPDH / G3PDH Protein (orb257495)

Description	Human GAPDH, His Tag is expressed from E. coli cells. It contains AA Gly 2 - Glu 335 (Accession # AA...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	GAPDH
Preservatives	50 mM Tris, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	N-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 36.9 kDa. The protein migrates as 37-38 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH01601
Purity	95%
Source	Human GAPDH, His Tag (orb257495) is expressed from E. coli cells. It contains AA Gly 2 - Glu 335 (Accession # AAH01601).
MW	36.9 kDa
Uniprot ID	P04406
NCBI	AAH01601
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human GAPDH protein