


## Product Datasheet

### Mouse GPA33 / A33 Protein (orb257518)

<b>Description</b>	Mouse GPA33, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 22 - Ile 235 (Ac...
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	GPA33
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-6×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.7 kDa. The protein migrates as 33-42 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_067623.1
<b>Purity</b>	95%
<b>Source</b>	Mouse GPA33, His Tag (orb257518) is expressed from human 293 cells (HEK293). It contains AA Leu 22 - Ile 235 (Accession # Q9JKA5-1).
<b>MW</b>	24.7 kDa
<b>Uniprot ID</b>	<a href="#">Q9JKA5</a>
<b>NCBI</b>	<a href="#">NP_067623.1</a>
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Mouse GPA33 protein