


Product Datasheet

Human HSP90AA1 / HSP86 Protein (orb257548)

Description	Human HSP90AA1, GST Tag is expressed from E. coli cells. It contains AA Glu 535 - Asp 732 (Accession...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	HSP90AA1
Preservatives	50 mM HEPES, 150 mM NaCl, pH7.0
Form/Appearance	Powder
Storage	-20°C
Tag	N-GST
Note	For research use only
Application notes	This protein carries a GST tag at the N-terminus. The protein has a calculated MW of 49.3 kDa. The protein migrates as 27 kDa,29 kDa,31-32 kDa and 47-49 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_005339
Purity	85%
Source	Human HSP90AA1, GST Tag (orb257548) is expressed from E.coli cells. It contains AA Glu 535 - Asp 732 (Accession # NP_005339).
MW	49.3 kDa
Uniprot ID	P07900
NCBI	NP_005339
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human HSP90AA1 prot...