


Product Datasheet

Human IgG2 Fc Protein (orb257567)

Description	Human IgG2 Fc, Tag Free is expressed from human 293 cells (HEK293). It contains AA Glu 99 - Lys 326 ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IgG2 Fc
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 25.7 kDa. The protein migrates as 33-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Glu 99 - Lys 326
Purity	90%
Source	Human IgG2 Fc, Tag Free (orb257567) is expressed from human 293 cells (HEK293). It contains AA Glu 99 - Lys 326 (Accession # P01859-1).
MW	25.7 kDa
Uniprot ID	P01859
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IgG Fc protei...