


Product Datasheet

Human IgG3 Fc Protein (orb257570)

Description	Human IgG3 Fc, Tag Free is expressed from human 293 cells (HEK293). It contains AA Glu 99 - Lys 377 ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IgG3 Fc
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 31.2 kDa. The protein migrates as 35-42 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Glu 99 - Lys 377
Purity	95%
Source	Human IgG3 Fc, Tag Free (orb257570) is expressed from human 293 cells (HEK293). It contains AA Glu 99 - Lys 377 (Accession # P01860).
MW	31.2 kDa
Uniprot ID	P01860
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IgG Fc protei...