


Product Datasheet

Mouse IgG2a Fc Protein (orb257575)

Description	Mouse IgG2a Fc, Tag Free is expressed from human 293 cells (HEK293). It contains AA Glu 98 - Lys 330...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IgG2a Fc
Preservatives	50 mM Tris, 100 mM Glycine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 26.4 kDa. The protein migrates as 30-35 kDa under reducing (R) condition, and 55-60 kDa under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	CAC20702
Purity	90%
Source	Mouse IgG2a Fc, Tag Free (orb257575) is expressed from human 293 cells (HEK293). It contains AA Glu 98 - Lys 330 (Accession # CAC20702).
MW	26.4 kDa
Uniprot ID	P01863
NCBI	CAC20702
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse IgG Fc protei...