



Product Datasheet

Human IL-4 Protein (orb257598)

Description	Human IL-4, premium grade is expressed from human 293 cells (HEK293). It contains AA His 25 - Ser 15...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	IL-4
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 15.0 kDa. The protein migrates as 19 kDa±3 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human IL-4, premium grade (orb257598) is expressed from human 293 cells (HEK293). It contains AA His 25 - Ser 153 (Accession # P05112-1).
MW	15.0 kDa
Uniprot ID	P05112
NCBI	AAH67514
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL4 protein



Immobilized Human IL-4, premium grade (C...