



## Product Datasheet

### Human IL-2 R beta / CD122 Protein (orb257617)

<b>Description</b>	Human IL-2 R beta, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ala 27 - Asp 23...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	IL-2 R beta
<b>Preservatives</b>	50 mM Tris, 100 mM Glycine, 150 mM NaCl, pH7.5
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.0 kDa. The protein migrates as 60-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Purity</b>	95%
<b>Source</b>	Human IL-2 R beta, Fc Tag (orb257617) is expressed from human 293 cells (HEK293). It contains AA Ala 27 - Asp 239 (Accession # NP_000869.1).
<b>MW</b>	51.0 kDa
<b>Uniprot ID</b>	<b>P14784</b>
<b>NCBI</b>	<b>NP_000869.1</b>
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human IL-2 R beta p...



Human IL-2 R beta, Fc Tag (Cat. No. orb2...