



## Product Datasheet

### Human IL-21 R Protein (orb257627)

<b>Description</b>	Human IL-21 R, His Tag is expressed from human 293 cells (HEK293). It contains AA Cys 20 - Pro 236 (...)
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	IL-21 R
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.9 kDa. The protein migrates as 37-47 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Purity</b>	90%
<b>Source</b>	Human IL-21 R, His Tag (orb257627) is expressed from human 293 cells (HEK293). It contains AA Cys 20 - Pro 236 (Accession # Q9HBE5-1).
<b>MW</b>	26.9 kDa
<b>Uniprot ID</b>	<a href="#">Q9HBE5</a>
<b>NCBI</b>	<a href="#">NP_068570.1</a>
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human IL-21 R prote...



Immobilized Biotinylated Human IL-21, Fc...