


Product Datasheet

Mouse gp130 / CD130 / IL-6 R beta Protein (orb257631)

Description	Mouse gp130, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Glu 617 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-6 R beta
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 93.3 kDa. The protein migrates as 120-130 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_034690
Purity	95%
Source	Mouse gp130, Fc Tag (orb257631) is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Glu 617 (Accession # NP_034690).
MW	93.3 kDa
Uniprot ID	Q00560
NCBI	NP_034690
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse gp130 protein