
Product Datasheet


Human Osteopontin (OPN) Protein (orb257724)

Description

Human Osteopontin, His Tag is expressed from



Reactivity	Human
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	Osteopontin
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 33.0 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	80%
Source	Human Osteopontin, His Tag (orb257724) is expressed from human 293 cells (HEK293). It contains AA Ile 17 - Asn 300 (Accession # AAH07016).
MW	33.0 kDa
Uniprot ID	P10451
NCBI	AAH07016
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Osteopontin p...